


Brain MRI processing for enhancing Alzheimer's diagnosis


Michela Proietti 1739846

Sofia Santilli 1813509



What is Alzheimer Disease (AD)?

- AD is a neurodegenerative disease that gradually affects neural cells, causing a progressive dementia
- It mainly affects the elderly population, even if initial signs of degeneration can be detected even in younger individuals
- The main indicators for the diagnosis of Alzheimer are tissue atrophy, precuneus and hippocampus atrophy



Magnetic resonance imaging (MRI)

- High spatial resolution
- Ability to contrast soft tissue
- Non-invasiveness
- Less health risks than PET or CT

Automatic analysis of brain MRI allows to have an **early diagnosis** of Alzheimer


- Artifacts
 - Respiratory and cardiac motion
 - Inappropriate spatial encoding
- Thermal noise
 - Hardware-induced errors
 - Body temperature

Accuracy limit in the current state-of-the-art for segmentation and classification methods

- Other noises
- Low SNR in some regions or tissues

**IMAGE
DENOISING**


K-space and noise distributions



State-of-the-art in MRI Denoising

- **Filtering, transformations, or statistical methods** (bilateral filtering, non-local means and BM3D currently are the most widely-used)
- **Convolutional NNs**: exploited for their small receptive field (SCNN, DnCNN)
- **Autoencoders**: may lose important information (as edge or fine structure)
- **U-net**: large receptive field by utilizing multi-scale features
- **GANs**: U-net often chosen as generator

Proposed approach



Our Novelty

□ Generator

- Denoising residual network
- Skip connections
- Batch normalization
- No pooling, stride or dilation

□ Discriminator

- Simple CNN
- It could be changed depending on the performance of the generator




Figure: Residual block

Validation and comparisons

1. Apply our denoising algorithm on different dataset for Alzheimer's classification
2. Use transfer learning in order to achieve a relatively fast but accurate classification
3. Compare the classification results obtained using our denoising algorithm to preprocess the images and those obtained without denoising

Expected results

- Increase in the classification accuracy by around 10% with respect to the classification done on a noisy dataset
- Not excessively long training time

References

- Daiki Tamada, *"Review: Noise and artifact reduction for MRI using deep learning"*. 28 February 2020, Department of Radiology, University of Yamanashi, Chuo, Yamanashi, 409-3898, Japan
- M. Tian and K. Song, *"Boosting Magnetic Resonance Image Denoising With Generative Adversarial Networks"*, in IEEE Access, vol. 9, pp. 62266-62275, 2021, doi: 10.1109/ACCESS.2021.3073944.
- Nagaraj Yamanakkanavar, Jae Young Choi and Bumshik Lee, *"Review: MRI Segmentation and Classification of Human Brain Using Deep Learning for Diagnosis of Alzheimer's Disease: A Survey"*, 7 June 2020, Department of Information and Communications Engineering, Chosun University, Gwangju 61452, Korea
- Haoyu Ren, Mostafa El-Khamy, Jungwon Lee, *"DN-ResNet: Efficient Deep Residual Network for Image Denoising"*, 16 October 2018
- Nikolas Adaloglou, *"Understanding the receptive field of deep convolutional networks"*, 2 July 2020
- Nikolas Adaloglou, *"Intuitive Explanation of Skip Connections in Deep Learning"*, 23 March 2020



Questions?